



MILLBROOK

# Investigation of a Degreasing Process using MiniSIMS / XPS

Surface residues are examined by a combination of analysis techniques to monitor the effectiveness of new aqueous degreasing agent.

As new environmental regulations come into force, changes to production processes are often required. One example is the move from organic solvents to aqueous (surfactant) based degreasing agents. Aqueous solutions pose fewer hazards during use and in subsequent disposal, but treatment parameters must be completely re-optimized for the new solutions. Surface analysis is required to monitor the effectiveness of the cleaning agent, and to detect any residues that could lead to later problems such as corrosion or poor adhesion of coatings.

Monitoring of the surface condition requires an analysis technique with high surface specificity and sensitivity. Both Secondary Ion Mass Spectrometry (SIMS) and X-ray Photoelectron Spectroscopy (XPS) are well-established methods, which give complementary information about the chemical composition of a surface. Each technique has its strengths; in particular, SIMS offers more information about molecular structure and a rapid imaging capability that is not available with lower cost XPS instrumentation. The extra information gained by visualization of the distribution of a contaminant can often give important clues about its origin.



This analysis examines a metallic component before and after cleaning with a proprietary aqueous degreasing agent. Residues are detected by XPS, but SIMS provides much more detailed information about their chemical identity and reveals the non-uniform distribution of some contaminants.


## Investigation of a Degreasing Process Using MiniSIMS / XPS

The component was first examined using an entry level (non-imaging) XPS instrument. As might be expected, the XPS spectrum of the uncleaned (blue) component shows a strong carbon peak. There are weak peaks corresponding to nickel and chromium, which are components of the metal alloy. This suggests the carbon layer is not continuous, and quantification of the XPS data shows carbon constitutes 76 atomic % of the surface layer.

The XPS spectrum after cleaning (red) shows more intense nickel and chromium peaks and a much diminished carbon signal. This indicates the cleaning process has been partially successful. However, sodium, phosphorus and silicon signals are clearly visible in the spectrum, indicating the presence of these elements at the surface. There is trace nitrogen signal due to cyano or nitride species.

The cleaned component was then examined using the MiniSIMS in surface spectrum mode (positive and negative SIMS). This confirmed the presence of a high level of sodium on the nickel / chromium surface. The chemical species containing silicon and phosphorus were revealed to be silicates and phosphates. The spectrum also confirms the nitrogen signal arises from a  $CN^-$  species at the surface. Other peaks in the spectra are due to the residual hydrocarbon species.






Chemical images of component

Above:  
Na<sup>+</sup> (Sodium)  
Cr<sup>+</sup> (Chromium)  
O<sup>-</sup> (Oxygen)

Below:  
CH<sup>-</sup> (Carbon)  
PO<sub>3</sub><sup>-</sup> (Phosphate)  
O<sup>-</sup> (Oxygen)

Image size  
2.5mm × 2.5mm



The major advantage of the MiniSIMS is that the distribution of these species can be mapped, an image being acquired in as little as 30 seconds. Most are seen to be uniform across the surface of the component. However the sodium is clearly present as non-uniform contamination. The shape suggests that this is a drying stain, proving that the sodium originates from the degreasing agent rather than arising from surface segregation in the metal itself.

Samples and XPS analysis courtesy of Dr. S Jenkins, LPD Labservices, LG Philips +44(0)1254 507377

MILLBROOK INSTRUMENTS LIMITED